

Surgical
DRAPES
GOWNS
PLASTIC
DISPOSABLE

tiomedikal.com

INDEX

TIOSET

QUICK CARE PACK
DRESSING PACK
QUICK CARE PACK
WOUND CARE PACK

GENERAL SURGERY

LAPAROTOMY PACK
LAPAROSCOPY PACK
UNIVERSAL PACK

CARDIOLOGY

VALVE PACK
ANGIOGRAPHY PACK
CARDIOVASCULAR PACK

BRAIN SURGERY

CRANIOTOMY PACK

OPHTHALMOLOGY

EYE PACK 1
EYE PACK 2

EAR - NOSE - THROAT

E.N.T PACK

GYNECOLOGY

DELIVERY PACK
HISTOSCOPY PACK
OOCYTE PICK-UP PACK
CESAREAN PACK
EMBRYO PACK

UROLOGY

PERCUTANEOUS PACK
CYCTOSCOPY PACK
T.U.R PACK

ORTHOPEDICES

HIP PACK
ARTROSCOPY PACK
EXTREMITY PACK

DRAPES

EYE DRAPE DOUBLE LCP
E.T DRAPE
T.U.R DRAPE
EYE DRAPE SINGLE LCP
EXREMITY DRAPE
CYSTOSCOPY DRAPE
"U" DRAPE
ANGIOGRAPHY D.
K.V.C DRAPE
CESAREAN DRAPE
LAZER EYE DRAPE
PERCUTANEOUS DRAPE
DRAPE W/ADHESIVE TAPE
CRANIOTOMY DRAPE
ADHESIVE TAPE
TRANSPARENT DRAPE
O.P.U TAPE
BABY DRAPE
PLAIN COVER
MAYO STAND DRAPE
LEG DRAPE
MEDICAL TOWEL
FOOT COVER
OP TAPE

TIOPLAS+

BOWLS

GALLIPOT 60ML
GALLIPOT 120ML
GALLIPOT 250ML
GALLIPOT 500ML
GALLIPOT 1000ML
GUIDEWIRE BOWL

TRAY

KIDNEY TRAY
TRAY 2500ML
TRAY 1000ML
THREE COMPARTENT TRAY 16*16CM

EQUIPMENTS

SPONGE FORCEPS 18/24CM
KOHER FORCEPS
SCALPEL
TWEEZERS
THREAD RETRIEVER
COIL REMOVAL
DILATOR 5/6-3/4MM
HYSTEROMETER
TENACULUM FORCEPS
TEALES VULSELUM FORCEPS
FORCEPS WITH METAL TIP

COVERS

ROUND TUBE COVERS
TUBE COVERS
ULTRASON PROPE COVER
MICROSCOPY COVER "ZEISS"
MICROSCOPY COVER "LECIA"
FLUOROSCOPY COVER
CAMERA COVER
C-ARM SCOPY COVER
RAW MATERIALS

SMS
SMMS
SPUNFLEX
DUBLEX
TRIPLEX
PP
"POLYPROPYLENE"
PE "POLIETILEN"
PU "POLIUURETAN"
GOWNS
REINFORCED SURGICAL GOWN
SURGICAL GOWN
FULL REINFORCED SURGICAL
GOWN PATIENT GOWN



SET

QUICK CARE PACK

CE ISO
13485 EN
13795



SET

• DRESSING PACK

537.01.R03.01



CONTENS	SIZE/CM	PCS
Three Compartment Tray	16*16	1
Non Woven Swabs	10*10	5
Non Woven Balls	3CM	4
Tweezers	11CM	3
Hand Towel	40*40	1
Drape	50*50	1
Pe Pag		1

• QUICK CARE PACK

538.01.R03.01



CONTENS	SIZE/CM	PCS
Solut on Cup	120ml	1
Solut on Cup	60ml	1
Non Woven Swabs	3CM	4
Fenestreted Drape	50*50	1
Tweezers Metal Tips	10,5CM	1
Tweezers Metal Tips	13,5CM	1
Tweezers Metal Tips	12CM	1
Pe Pag		1
Hand Towel	40*40	1
Drape	50*50	1

• WOUND CARE PACK

539.01.R03.01



CONTENS	SIZE/CM	PCS
Kidney Tray	700ml	1
Sponge Forceps	18cm	1
Nitrile Gloves	M	2
Gauze Swabs	10*10	3
Fenestreted Drape	50*50	1
Hand Towel	40*40	1
Solut on Cup	250ml	1
Gauze Ball	3CM	4
Plaster	9*10	1

GENERAL SURGERY

CE ISO
13485 EN
13795


TiO
SET

• UNIVERSAL SET

510.OX.R03.01

Intended use and place of use:


The product is intended for general purposes. All surgical procedures can be performed with this set. The set includes covers that cover the top and sides of the patient. When the anesthesia drape, side band drape and foot drape are joined together, a drape with the surgical area in the center is created on the patient.



CONTENS	SIZE/CM	PCS
Anesthesia Drape	150*240	1
Leg Drape	150*180	1
Drape W/Adhesive Tap	75*90	2
Mayo Stand Drape	75*145	1
Back Table Cover	150*200	1
Reinforced Surgical Gown	L	3
Medical Towel	40*40	3
Op-Tape	10*30	2

50*80*50	40*60*40
PIECES IN BOX	PIECES IN BOX
20	10
HEAVY	HEAVY
1115gr	1115gr

FABRIC OPTION



FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Highly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




GENERAL SURGERY

LAPAROTOMY SET

511.OX.R03.01

Intended use and place of use:


This product is intended for laparotomy surgery. There is a window on the drape, which is covered with a polyurethane incision film. There are also PE fluid collection bags around the window. These products are medical devices that may only be used once for a single patient



CONTENS	SIZE/CM	PCS
Laparotomy Drape	250*300	1
Back Table Cover	150*200	1
Reinforced Surgical Gown	L	3
Medical Towel	40*40	3
Op Tape	10*30	3

50*80*50	40*60*40
PIECES IN BOX	PIECES IN BOX
10	5
HEAVY	HEAVY
1200gr	1200gr

FABRIC OPTION



Medical
nonwoven
fabric

Medical blue

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




ALL SIZES ARE PRODUCED & CUSTOMIZATION

• LAPAROSCOPY SET

512.OX.R03.01

Intended use and place of use:


This product is intended for laparoscopic surgery. On the drape there is a window covered with polyurethane incision film. There are also PE fluid collection bags around the window.



CONTENS	SIZE/CM	PCS
Laparotomy Drape	250*300	1
Back Table Cover	150*200	1
Mayo Stand Drape	75*145	1
Reinforced Surgical Gown	L	3
Medical Towel	40*40	3
Op Tape	10*30	3

50*80*50	40*60*40
PIECES IN BOX	PIECES IN BOX
10	5
HEAVY	HEAVY
1200gr	1200gr

FABRIC OPTION



FABRIC

Medical
nonwoven

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Highly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• ANGIOGRAPHY SET

553.OX.R03.01

Intended use and place of use:


It is used for angiography operations. There are 2 holes which are 70 cm from the head. The holes are covered with a polyurethane cutting film. The fluoroscopy sheath included in the set is used to protect the existing device from contamination during surgery.



CONTENS	SIZE/CM	PCS
Angiography Drape	150*300	1
Back Table Cover	100*150	1
Gavze Swabs	7.5*7.5	10
Sponge Stick	15	1
Fluoroscopy Cover	86	1
Solution Cup	250cc	2
Surgical Gown	L	2
Medical Towel	40*40	2

50*80*50	40*60*40
PIECES IN BOX	PIECES IN BOX
20	10
HEAVY	HEAVY
778gr	778gr

FABRIC OPTION



FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• BY-PASS (CARDIOVASCULAR) DRAPE SET

551.OX.R03.01

Intended use and place of use:

The product is used during open heart surgery. The cardiovascular dressing is placed over the patient, following the directional instructions on the dressing. The patient has a perforated window in the chest area, which is covered with an incision drape. The kit includes a special U-slot Cover that covers the perineal area. There are 2 foot supports that are adapted to the patient's feet. There are waste bags for different waste.



CONTENS	SIZE/CM	PCS
By-Pass Drape	240*340	1
"U"Drape w/Perine Drape	150*240	1
Back Table Cover	150*200	1
Mayo Stand Drape	74*145	1
Reinforced Surgical Gown	L	3
Defibulator Bag	30*35	2
Op-Tape	10*30	2
Foot Cover	standart	2
Medikal Towel	40*40	3

50*80*50

40*60*40

PIECES IN BOX

PIECES IN BOX

10

5

HEAVY
2050gr

HEAVY
2050gr

FABRIC OPTION



SMS



DUBLEX



TRIPLEX



PP

FABRIC

COLOR

WEIGHT

PROPERTIES

Medical
nonwoven
fabric

Medical blue

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
TRIPLEX 64-70 GR
SPUNFLEX 50-54

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• VALVE-SET

552.OX.R03.01

Intended use and place of use:

The product is used during open heart surgery. The cardiovascular dressing is placed over the patient, following the directional instructions on the dressing. There is a perforated window in the patient's chest area that is covered with an incision drape. There are 2 foot covers for the patient's feet. There are waste bags for various wastes.



CONTENS	SIZE/CM	PCS
K.V.C Drape	240*340	1
Back Table Cover	150*200	1
Mayo Stand Drape	75*145	1
Foot Cover	36*45	2
Defibillator Bag	30*35	2
Op-Tape	10*30	2
Reinforced Surgical Gown	L	3
Medical Towel	40*40	1

50*80*50

40*60*40

PIECES IN BOX

10


HEAVY
1540gr

PIECES IN BOX

5

HEAVY
1540gr

FABRIC OPTION



FABRIC

Medical
nonwoven

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
TRIPLEX 64-70 GR
SPUNFLEX 50-54

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic



BRAIN SURGERY

ISO
13485 EN
13795




• CRANIOTOMY SET

571.OX.R03.01

Intended use and place of use:

The craniotomy cover set is used for neurosurgical procedures. There is a fluid collection bag on the craniotomy drape and the drape completely covers the patient. The surgical area is covered with a polyurethane incision drape. The fluid collection bag is located under the incision area. Side strip covers are placed on the sides of the patient.



CONTENS	SIZE/CM	PCS
Craniotomy Drape	200*350	1
"U" Drape	150*250	1
Drape W/Adhesive Tape	75*90	4
Mayo Stand Drape	75*145	1
Back Table Cover	150*200	1
Reinforced Surgical Gown	L	3
Medical Towel	40*40	3

50*80*50

40*60*40

PIECES IN BOX

PIECES IN BOX


15

7

HEAVY
1647gr

HEAVY
1647gr

FABRIC OPTION



FABRIC

COLOR

WEIGHT

PROPERTIES

Medical
nonwoven
fabric

Medical blue

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
TRIPLEX 64-70 GR
SPUNFLEX 50-54

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• EYE PACK 1

591.OX.R03.01

Purpose and location of use:

The product is used for eye surgeries. There is 1 fluid collection bag in the eye patch. The eye cover can only be used for applications for one eye. The flat cloth is used to cover the rest of the patient.



CONTENS	SIZE/CM	PCS
Eye Drape Single LCP	100*150	1
Plain Drape	150*180	1
Back Table Cover	100*150	1
Surgical Gown	L	2

50*80*50

40*60*40

PIECES IN BOX

40

HEAVY
500gr

PIECES IN BOX

20

HEAVY
500gr

FABRIC OPTION



SMS



DUBLEX



TRIPLEX



PP

FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic





• EYE PACK 2

592.OX.R03.01

Intended use and place of use:

The product is used for eye surgery. There are 2 fluid collection bags (pouches) in the eye cover. The eye cover can be used for both eyes. The flat cloth is used to cover the rest of the patient.



CONTENS	SIZE/CM	PCS
Eye Drape Double LCP	100*150	1
Plain Drape	150*180	1
Back Table Cover	100*150	1
Surgical Gown	L	2
Medical Towel	40*40	2

50*80*50

40*60*40

PIECES IN BOX

PIECES IN BOX

35

17

HEAVY
600gr

HEAVY
600gr

FABRIC OPTION



SMS



DUBLEX



TRIPLEX



PP

FABRIC

COLOR

WEIGHT

PROPERTIES

Medical
nonwoven
fabric

Medical blue

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic



EAR, NOSE & THROAT

ISO
13485 EN
13795




• ENT SET

560.OX.R03.01

Purpose and location of use:

The product is used in operations in the throat, nose, and ear area. The ENT cloth has a U-shaped slit and is surrounded by adhesive tape. The flat cloth is used to cover the rest of the patient.



CONTENS

E.N.T Drape
Plain Drape
Back Table Cover
Surgical Gown
Medical Towel

SIZE/CM

150*250
100*100
100*150
L
40*40

PCS

1
1
1
1
2

50*80*50

40*60*40

PIECES IN BOX

20


HEAVY
824gr

PIECES IN BOX

10

HEAVY
824gr

FABRIC OPTION



FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• CESAREAN SET

560.OX.R03.01

Purpose and place of use:


The product is manufactured for cesarean section surgeries. The drape has a 25 x 25 cm window and is covered with a polyurethane cutting film. Around the window, there is a fluid collection bag. On the fluid collection bag (pouch), there is an aspiration plug. Baby blanket, for wrapping the newborn.



CONTENS	SIZE/CM	PCS
Cesarean Drape	200*300	1
Back Table Cover	150*200	1
Baby Drape	60*80	1
Reinforced Surgical Gown	L	3
Medical Towel	40*40	2

50*80*50	40*60*40
PIECES IN BOX	PIECES IN BOX
20	10
HEAVY	HEAVY
1127gr	1127gr

FABRIC OPTION



FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• DELIVERY-SET

541.OX.R03.01

Purpose and location of use:


The product is used in delivery operations. The layers of the bundled cloth are unfolded on the patient. There are 2 leg covers that are placed separately over the legs. The flat blanket is placed under the patient. The baby blanket is used to wrap the newborn. The midwives wear surgical gowns.



CONTENS	SIZE/CM	PCS
Drape W/Adhesive Tape	75*150	1
Plain Drape	75*100	1
Baby Drape	60*80	1
Back Table Cover	100*150	1
Surgical Gown	L	2
Medical Towel	40*40	2
Leg Cover	60*120	2

50*80*50	40*60*40
PIECES IN BOX	PIECES IN BOX
25	13
HEAVY	HEAVY
815gr	815gr

FABRIC OPTION



FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• OOCYTE PICK-UP PACK

543.OX.S03.01

Purpose and location of use:

The product is used for IVF operations. In the lithotomy position, 2 leg cuffs are separately attached to the patient's legs. The set includes an OPU cover and a trapezoidal hole surrounded by adhesive tape. The adhesive tape-covered hole part is fixed in the perineal area and the operation is performed.



CONTENS

O.P.U Drape

Plain Drape

Transparent Drape

Back Table Cover

Leg Cover

Surgical Gown

Medical Towel

SIZE/CM

130*150

75*100

60*60

100*150

60*120

L

40*40

PCS

1

1

1

1

2

2

2

50*80*50

40*60*40

PIECES IN BOX

25

HEAVY
815gr

PIECES IN BOX

13

HEAVY
815gr

FABRIC OPTION



SMS



DUBLEX



TRIPLEX



PP

FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• EMBRYO-SET

546.0X.000.01

Purpose and place of use:

It is used for embryo transfer operations. In the middle of the E.T. cover there is a hole surrounded by adhesive tape. It is placed under the patient with the flat cloth.



CONTENS

E.T Drape
Plain Drape
Gauzes

SIZE/CM

80*90
75*90
7.5*7.5

PCS

1
1
10

50*80*50

40*60*40

PIECES IN BOX

150
HEAVY
103gr

PIECES IN BOX

73
HEAVY
103gr

FABRIC OPTION



SMS



DUBLEX



TRIPLEX



PP

FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

PROPERTIES


Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic



• T.U.R PACK

521.03.R02.01

The product is manufactured for urological procedures and used in male patients.



CONTENS	SIZE/CM	PCS
T.U.R Drape	175*240	1
Plain Drape	100*150	1
Back Table Cover	150*200	1
Reinforced Surgical Gown L		2
Medical Towel	40*40	2

50*80*50

40*60*40

PIECES IN BOX

20

HEAVY
910gr

PIECES IN BOX

10

HEAVY
910gr

FABRIC OPTION



SMS



DUBLEX



TRIPLEX



PP

FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• PERCUTANEOUS PACK

522.01.R02.01

Purpose and location of use:

The product is used in kidney surgeries. In the center of the drape is a window measuring 25 x 25 cm, which is covered with a polyurethane incision film. Under the incision area is a fluid collection bag. The percutaneous cover set includes 1 camera cover, which keeps the camera sterile.



CONTENS	SIZE/CM	PCS
Percutan Drape	150*300	1
Back Table Cover	150*200	1
Op-Tape	10*30	1
Scopy Cover	60*60	1
Camera Cover	14*250	1
Reinforced Surgical Gown	L	2
Medical Towel	40*40	2

50*80*50

40*60*40

PIECES IN BOX

PIECES IN BOX


20

10

HEAVY
916gr

HEAVY
916gr

FABRIC OPTION



FABRIC

COLOR

WEIGHT

PROPERTIES

Medical
nonwoven
fabric

Medical blue

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• CYSTOSCOPY PACK

523.03.R02.01

Purpose and location of use:

The cover set for cystoscopy is used in gynecological procedures on female patients. The patient's legs are covered with the help of leg covers. The cystoscopy cover set is placed on the perineal area of the patient and the patient is prepared for the procedure. The flat cover is placed under the patient.



CONTENS	SIZE/CM	PCS
Cystoscopy Drape	75*150	1
Plain Drape	100*120	1
Back Table Cover	150*200	1
Reinforced Surgical Gown L		2
Leg Cover	60*120	1
Medical Towel	40*40	2

50*80*50

40*60*40

PIECES IN BOX

PIECES IN BOX


30

15

HEAVY
734gr

HEAVY
734gr

FABRIC OPTION



SMS

DUBLEX

TRIPLEX

PP

FABRIC

COLOR

WEIGHT

PROPERTIES

Medical
nonwoven
fabric

Medical blue

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• HIP PACK

534.01.R03.01

Purpose and location of use:

The product is used for hip surgical procedures. It consists of a U-slot cover. The U-slot cover is placed between the legs, exposing the patient's hip, and secured to the surgical area using adhesive tape. Then the anesthesia drape is pulled over and the patient is prepared. Additionally, flat drapes are used to cover the patient.



CONTENS	SIZE/CM	PCS
"U" Drape	200*300	1
Anesthesia Drape	150*240	1
Plain Drape 1	150*200	1
Plain Drape 2	75*90	1
Back Table Cover	150*200	1
Mayo Stand Drape	75*145	1
Op-Tape	10*50	2
Medical Towel	40*40	3
Reinforced Surgical Gown L		3

50*80*50

40*60*40

PIECES IN BOX

PIECES IN BOX


10

5

HEAVY
1404gr

HEAVY
1404gr

FABRIC OPTION



SMS

DUBLEX

TRIPLEX

PP

FABRIC

COLOR

WEIGHT

PROPERTIES

Medical
nonwoven
fabric

Medical blue

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




■ ARTHROSCOPY-SET

531.OX.R03.01

Use and location of use:

The product is used in knee surgeries involving exposure to liquids. The bag-like arthroscopy drape is used to cover the patient. The drape has 2 elastic parts. This part is used to comfortably guide the patient's knee into the bag part. The surgery is performed and up to 20 liters of fluid are retained by the special structure of the bag. Excess fluid is suctioned out of the system using the suction plug. The leg cuff is used to wrap around the patient's leg. The set includes a camera cover that keeps the camera sterile.



CONTENS	SIZE/CM	PCS
Arthroscopy Drape	200*300	1
Back Table Cover	150*200	1
Camera Cover	14*250	1
Reinforced Surgical Gown	L	3
Leg Cover	36*45	1
Op-Tape	10*30	1
Medical Towel	40*40	2

50*80*50

40*60*40

PIECES IN BOX

20


HEAVY
1184gr

PIECES IN BOX

10

HEAVY
1184gr

FABRIC OPTION



FABRIC

COLOR

WEIGHT

PROPERTIES

Medical
nonwoven
fabric

Medical blue

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic




• EXTREMITY PACK

536.02.R03.01

Purpose and location of use:

The product is used in knee surgeries. The limb cover is used to cover the patient. The cover is equipped with an elastic part with a hole. The patient's leg is passed through this hole before the surgery.



CONTENS	SIZE/CM	PCS
Exremity Drape	200*300	1
Back Table Cover	150*200	1
Op-Tape	10*30	2
Mayo Stand Drape	75*145	1
Skin Marker&Ruler		1
Reinforced Surgical Gown	L	3
Medical Towel	40*40	2

50*80*50

PIECES IN BOX

20

HEAVY
1096gr

40*60*40

PIECES IN BOX

10

HEAVY
1096gr

FABRIC OPTION



SMS

DUBLEX

TRIPLEX

PP

FABRIC

Medical
nonwoven
fabric

COLOR

Medical blue

WEIGHT

SMS/SMMS 35-43GR
DUBLEX 52-56 GR
SPUNFLEX 50-54 GR

PROPERTIES

Antistatic
Opaque
Soft touch
Hydrophobic
Hyghly Breathable
Metal-free
Impervious to microbial
migration
OEKO-TEX certified
ISO 9001
Antiallergic





DRAPES



SMS

It consists of 3 separate shafts, which are connected thermally or with adhesive.



SMMS

It consists of 4 separate shafts, which are either thermally or adhesively connected.



SPUNFLEX

The product is made by laminating the chipboard with the foil. Since it is not very absorbent, it is only used for the back wall.



DUBLEX

It consists of a polyethylene film and an absorbent material (chemical band). The absorbent part has a higher absorbency than the nonwoven fabric.



TRIPLEX

It consists of 3 materials. It is formed by the lamination of absorbent material, polyethylene film, and spunlace fabric.



PP "Polypropylen"



PE "Polyethylen"


DRAPES



ISO
13485 EN
13795




546.03.001.01



E.T. Cover
SIZE/CM

A	B
90	80
100	100


521.03.R02.01



T.U.R COVER
SIZE/CM

A	B
240	175
260	200


441.03.001.01



EYE COVER
CLOTH
SIMPLE LCP
SIZE/CM

A	B
100	130
100	150
150	150


523.03.R.02.01



CYSTOSCOPY-
DRAPE
SIZE/CM

A	B
150	75

536.02.R03.01



EXTREMITIES
COVER
SIZE/CM

A	B
300	200
320	255
330	230


571.02.000.01



CRANIOTOMY-
COVER
SIZE/CM


A	B
350	200
290	230

ALL SIZES ARE PRODUCED & CUSTOMIZATION



DRAPES


534.01.000.01



"U" COVER
SIZE/CM

A	B
300	200
250	150


553.01.001.01



ANGIOGRAPHY-
COVER
SIZE/CM 150*300

A	B
300	150
330	220

551.02.R03.01




K.V.C DRAPE

SIZE/CM

A	B	C
340	240	300
350	220	300


542.01.R03.01



CESAREAN -
ABDOMINAL COVER
SIZE/CM

A	B
300	200

441.03.001.01




EYE LASER

SIZE/CM

A	B
60	40
41	41

522.01.R02.01



PERCUTAN COVERING
SIZE/CM

A	B
300	150
150	200

ALL SIZES ARE PRODUCED & CUSTOMIZATION




DRAPES



ISO
13485 EN
13795



441.03.001.01




DOUBLE LCP
EYE COVER

SIZE/CM

A	A	B
100		130
100		150
150		150

543.03.S02.01



O.P.U COVER

SIZE/CM 130*150

A	B
150	130
200	240
150	200

ASSOCIATION SIMPLE



COVER CLOTH


SIZE/CM 60*60



BABY TOWEL

SIZE/CM 60*80

445.03.006.01



DRESS TOWEL
WITH ADHESIVE
STRIPS

SIZE/CM

A	B	C
50	50	40
60	60	50
75	75	65
100	100	90



PLAIN DRAPE
SIZE/CM

A	B
50	50
60	60
75	75
100	100



SMS



DUBLEX



TRIPLEX



PP

ALL SIZES ARE PRODUCED & CUSTOMIZATION

ASSOCIATION SIMPLE

CE ISO
13485 EN
13795




444.03.007.01



PLAIN DRAPE
SIZE/CM

A	B
150	100
150	150
200	150
250	200
190	150

448.03.000.01



MAYO TABLE DRAPE


SIZE/CM

A	B
145	75
145	79
145	80



LEG COVER
SIZE/CM

A	B
120	60




HAND TOWEL
SIZE/CM

A	B
20	40
40	40
30	40



FOOT COVER

SIZE/CM 30*35



OP-TAPE
SIZE/CM

A	B
50	10
30	10
15	10



SMS



DUBLEX



TRIPLEX



PP



TiOPLAS+

BOWLS

CE ISO
13485 EN
13795



GALLIPOT
622.01.000.01
VOLUME
60ml



GALLIPOT
622.01.001.01
VOLUME
120ml



GALLIPOT
622.01.002.01
VOLUME
250ml



GALLIPOT
622.01.003.01
VOLUME
500ml



GALLIPOT
622.01.004.01
VOLUME
1000ml



**GUIDEWIRE
BOWL**
622.01.004.01
VOLUME
2500ml

ALL COLOURS ARE PRODUCED



SPONGE FORCEPS

640.01.000.01

VOLUME

24cm

18cm



TENACULUM

FORCEPS

640.01.006.01

VOLUME

26cm



TWEEZERS

640.01.002.01

VOLUME

13cm



KOHER FORCEPS

640.01.003.01

VOLUME

14cm




COIL REMOVAL

640.01.005.01

VOLUME

25cm



**TEALES
VULSELUM
FORCEPS**

660.01.000.01

VOLUME

27cm

EQUIPMENTS

CE ISO
13485 EN
13795



**FORCEPS WITH
METAL TIP
640.01.007.01**

VOLUME
12,5cm



**FORCEPS WITH
METAL TIP
640.01.008.01**

VOLUME
10,5cm



**FORCEPS WITH
METAL TIP
640.01.009.01**

VOLUME
13,6cm



**SCALPEL
640.01.004.01**

VOLUME
Var ous Numbers



**THREAD RETRIEVER
650.01.002.01**

VOLUME
25cm



**DILATOR
650.01.000.01**

VOLUME
5/6mm
3/4mm




**HYSTEROMETER
650.01.003.01**

VOLUME
25cm


ALL COLOURS ARE PRODUCED

TUBE COVERS


Round Tube Covers 219.10.001.01




Tube Covers 219.10.003.01



C-Arm Cover W th Cl ps 219.10.006.01



ULTRASOUND PROBE COVERS

PRODUCT/CODE	PURPOSE OF USAGE	TYPE	WIDTH	LENGTH
 219.10.007.01	General Purpose	PP/PU Telescopically Folded Latex Free Titled and straight (Optionally, in kit form with Elastic Band and 20ml Ultrasound gel)	100mm/3,93" 120mm/4,72" 140mm/5,51" 150mm/5,90" 160mm/6,30" 180mm/7,08"	1000mm/39,37" 1220mm/48,03" 1440mm/55,11" 2440mm/96,06"
 219.10.008.01	Obstetrics Gynecology	PU Flat Folded Latex Free Titled and straight (Optionally, in kit form with Elastic Band and 20ml Ultrasound gel)	52mm*127mm 2,04"*5,00"	610mm/24,01"
 219.10.009.01	Surgical	PP/PU Telescopically Folded Latex Free Titled (Optionally, in kit form with Elastic Band and 20ml Ultrasound gel)	130mm-10mm/5"-0,39" 130mm-20mm/5"-0,79" 130mm-40mm/5"-1,57" 150mm-90mm/5,90"-3,54" 140mm/5,51"	2440mm/96,06"
 219.10.010.01	Trans Esophageal	PU Soft Probe Cover (On applicator), Bite Block, Sterile Gel in Syringe and tapes	50mm-30mm 1,96"-1,18"	1220mm/48,03"
 219.10.011.01	Other	PP/PU Telescopically Folded Latex Free Titled and Straight (Optionally, in kit form with Elastic Band and 20ml Ultrasound gel)	80mm/3,15" 100mm/3,93" 140mm/5,51" 180mm/7,08"	300mm/11,81" 610mm/24,07" 1200mm/47,24"

ALL SIZES ARE PRODUCED & CUSTOMISATION

TEE/TOE PROPE COVER KIT

219.10.012.01



INDEX

- Tee/Toe Prope Cover W th
Ultrasound Gel Syringe
- K dney Tray
- Endoscopy Mouthp ece
- Apron
- Wrap
- T e 2/ pcs

Specially designed for Trans Esophageal or Trans Esophageal Echocardiography (TOE / TEE) transducers for adults and neonatal patients. As TIO MEDICAL, we have prepared an easy-to-use, complete k t w th probe protector, ultrasound gel, applicator, apron and cup. TIO MEDICAL offers a complete range of Probe protection kits to protect pat ents and healthcare workers against cross contamination for TOE / TEE procedures.

All TOE / TEE ultrasound transducer covers are made of latex-free PU material.

All of our TOE/TEE ultrasound transducer sleeves will f t every commeri ally available TEE transducer, adult, pediatric and Micromulti .



MICROSCOPY COVER

"ZEISS"

221.10.000.01

Mater al: PE



PE



MICROSCOPY COVER

"LECIA"

222.10.000.02

Mater al: PE



PE



SCOPY COVER

219.10.000.00

Mater al: PE



PE

SIZE/CM:
60*60cm
80*80cm
90*90cm
130*130cm
100*100cm
120*120cm
120*100cm
80*120cm



ULTRASOUND
PROBE COVER

219.10.010.01

Mater al: PE



PE

SIZE/CM:
15*100
15*150
10*61



CAMERA COVER

211.10.000.01

Mater al: PE



PE

SIZE/CM:
250*14cm



C-ARM SCOPY
COVER

219.10.006.01

Mater al: PE



PE

SIZE/CM:
40*120cm

FULL REINFORCED SURGICAL GOWN



Outer View



Inter or View



Small	331.03.F01.01
Medium	332.03.F01.01
Large	333.03.F01.01
X-Large	334.03.F01.01
Xx-Large	335.03.F01.01

Material: SMS



SMS



REINFORCED SURGICAL GOWN



Outer View



Inner or View



Small	311.03.R01.01
Medium	312.03.R01.01
Large	313.03.R01.01
X-Large	314.03.R01.01
Xx-Large	315.03.F01.01

Material: SMS



SMS



STANDART GOWN



Small	321.03.S01.01
Medium	322.03.S01.01
Large	323.03.S01.01
X-Large	324.03.S01.01

Material: SMS



SMS



PATIENT GOWN



Small	351.06.P01.02
Medium	352.06.P01.02
Large	353.06.P01.02
X-Large	354.06.P01.02

Material: SMS



SS



Together **STRONGER**



TİO MEDICAL PRODUCTS INDUSTRY AND TRADE LTD.COMP.
TİO MEDİKAL TIBBİ ÜRÜNLER TUR. İNŞ. SAN. TİC. LTD. ŞTİ.
2/20 sok. No.:53 (Begos 3, Kuzey Giriş, 35400 Buca Osb/Buca/İzmir
Phone : +90 232 484 33 11 - 431 09 01 Fax : +90 232 484 33 22
nfo@t.omed.kal.com - export1@t.omed.kal.com
www.t.omed.kal.com



tiomedikal_



ISO
13485 EN
13795

